

Resource Summary Report

Generated by [dkNET](#) on Aug 16, 2026

Leica: ThermoBrite Slide Denaturation/Hybridization System

RRID:SCR_028505

Type: Tool

Proper Citation

Leica: ThermoBrite Slide Denaturation/Hybridization System (RRID:SCR_028505)

Resource Information

URL: <https://www.leicabiosystems.com/ihc-ish/fish-molecular-solutions/thermobrite/>

Proper Citation: Leica: ThermoBrite Slide Denaturation/Hybridization System (RRID:SCR_028505)

Description: Programmable, automated laboratory instrument used for slide based FISH (Fluorescence In Situ Hybridization) and ISH procedures. It is specialized, microprocessor controlled hot plate designed specifically for the denaturation and hybridization of DNA on glass microscope slides.

Synonyms: , Leica: ThermoBrite System, Leica ThermoBrite System

Resource Type: instrument resource

Keywords: denaturation and hybridization of DNA, glass microscope slides, slide based FISH, Fluorescence In Situ Hybridization, ISH procedures, microprocessor controlled hot plate,

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR_028505.pdf

Resource Name: Leica: ThermoBrite Slide Denaturation/Hybridization System

Resource ID: SCR_028505

Record Creation Time: 20260605T195314+0000

Record Last Update: 20260815T183315+0000

Ratings and Alerts

No rating or validation information has been found for Leica: ThermoBrite Slide Denaturation/Hybridization System.

No alerts have been found for Leica: ThermoBrite Slide Denaturation/Hybridization System.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.