

Resource Summary Report

Generated by [dkNET](#) on Sep 7, 2026

Ethanol precipitation of RNA from small or large volumes

DOI:10.17504/protocols.io.36fgrbn

Type: Protocol

Proper Citation

Stephen Floor 2019. Ethanol precipitation of RNA from small or large volumes.
protocols.io dx.doi.org/10.17504/protocols.io.36fgrbn

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.36fgrbn>

Authors: Stephen Floor

Group: Stephen Floor Lab

Summary: How to precipitate RNA with ethanol and resuspend it! Two separate protocols - one for small and one for large volumes. The procedure for large volume precipitation is similar, but begins with an initial precipitation and transfers to a small tube for the final steps to avoid loss in the large tube. Strongly recommended to move to a smaller tube. Recommended reading: Walker & Lorsch RNA purification--precipitation methods.

Affiliations: UCSF

Version: 1

Publication Date: 2019

Proper Citation: Stephen Floor 2019. Ethanol precipitation of RNA from small or large volumes. protocols.io dx.doi.org/10.17504/protocols.io.36fgrbn

Record Creation Time: 20210329T220529+0000

Record Last Update: 20210329T220937+0000

Data and Source Information

Source: