

Resource Summary Report

Generated by [dkNET](#) on Sep 14, 2026

MojoSort™ Streptavidin Nanobeads Column Protocol - Positive Selection

DOI:10.17504/protocols.io.692hh8e

Type: Protocol

Proper Citation

Kelsey Miller, Sam Li 2019. MojoSort™ Streptavidin Nanobeads Column Protocol - Positive Selection. protocols.io dx.doi.org/10.17504/protocols.io.692hh8e

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.692hh8e>

Authors: Kelsey Miller, Sam Li

Group: BioLegend

Summary: BioLegend MojoSort™ nanobeads work in commonly used separation columns, based on our internal research as well as validation by external testing by academic labs. This simple protocol consists of following the MojoSort™ protocol to label the cells with pre-diluted MojoSort™ reagents and using the columns as indicated by the manufacturer. Note: Due to the properties of our beads, it may be possible to use far fewer beads and less antibody than with other commercial suppliers. We recommend a titration to find the best dilution factor. However, as a general rule, dilutions ranging from 1:3 to 1:20 for the Nanobeads can be used. Please contact BioLegend Technical Service (tech@biolegend.com) if further assistance is needed.

Affiliations: BioLegend, BioLegend

External URL: <https://www.biolegend.com/protocols/mojosort-streptavidin-nanobeads-column-protocol-positive-selection/4773/>

Version: 2

Publication Date: 2019

Proper Citation: Kelsey Miller, Sam Li 2019. MojoSort™ Streptavidin Nanobeads Column Protocol - Positive Selection. protocols.io dx.doi.org/10.17504/protocols.io.692hh8e

Record Creation Time: 20210329T220521+0000

Record Last Update: 20210329T220849+0000

Data and Source Information

Source: