

Resource Summary Report

Generated by [dkNET](#) on Sep 14, 2026

MojoSort™ Streptavidin Nanobeads Protocol

DOI:10.17504/protocols.io.e2rbgd6

Type: Protocol

Proper Citation

Kelsey Miller 2016. MojoSort™ Streptavidin Nanobeads Protocol. protocols.io
[dx.doi.org/10.17504/protocols.io.e2rbgd6](https://doi.org/10.17504/protocols.io.e2rbgd6)

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.e2rbgd6>

Authors: Kelsey Miller

Group: BioLegend

Summary: The cells targeted by the Streptavidin Nanobeads are either selected or depleted by incubating your sample with the magnetic particles after incubating with a biotin-conjugated antibody or antibody cocktail. The magnetically labeled fraction is retained by the use of a magnetic separator. After collection of the targeted cells, downstream applications include functional assays, gene expression, phenotypic characterization, etc.

Affiliations: BioLegend

External URL:

http://www.biolegend.com/media_assets/support_protocol/MojoSort_Streptavidin_Nanobeads_Protocol

Version: 1

Publication Date: 2016

Proper Citation: Kelsey Miller 2016. MojoSort™ Streptavidin Nanobeads Protocol. protocols.io [dx.doi.org/10.17504/protocols.io.e2rbgd6](https://doi.org/10.17504/protocols.io.e2rbgd6)

Record Creation Time: 20210329T220521+0000

Record Last Update: 20210329T220848+0000

Data and Source Information

Source: