

Resource Summary Report

Generated by [dkNET](#) on Sep 16, 2026

Flex-T™ HLA Class I ELISA Protocol

DOI:10.17504/protocols.io.e7ebhje

Type: Protocol

Proper Citation

Kelsey Miller 2016. Flex-T™ HLA Class I ELISA Protocol. protocols.io
[dx.doi.org/10.17504/protocols.io.e7ebhje](https://doi.org/10.17504/protocols.io.e7ebhje)

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.e7ebhje>

Authors: Kelsey Miller

Group: BioLegend

Summary: The HLA class I ELISA is an enzyme immunoassay based on the detection of ?2-microglobulin subunit of HLA class I complexes, after capturing the complex through the conjugated biotin. To this end, biotinylated HLA class I complex is first captured in streptavidin coated microtiter wells. Subsequently, HRP-conjugated anti-human ?2-microglobulin is added to detect intact HLA class I complexes. Only intact HLA class I complexes are recognized. Peptides with high affinity binding will be clearly detected by this ELISA technique, while peptides with a moderate to low binding affinity for HLA class I provide a moderate to non-detectable signal. This protocol is designed to evaluate the efficiency of peptide exchange when using the Flex-T™ system.

Affiliations: BioLegend

External URL:

http://www.biolegend.com/media_assets/support_protocol/Protocol%20for%20HLA%20class%20I%20ELISA%20Protocol.pdf

Version: 1

Publication Date: 2016

Proper Citation: Kelsey Miller 2016. Flex-T™ HLA Class I ELISA Protocol. protocols.io
[dx.doi.org/10.17504/protocols.io.e7ebhje](https://doi.org/10.17504/protocols.io.e7ebhje)

Record Creation Time: 20210329T220516+0000

Record Last Update: 20210329T220830+0000

Data and Source Information

Source: