

Resource Summary Report

Generated by [dkNET](#) on Sep 8, 2026

ChroPlate - ProteinG

DOI:10.17504/protocols.io.uyjexun

Type: Protocol

Proper Citation

Alexandra Ehl, David Frommholz, Nadine Stefanczyk 2018. ChroPlate - ProteinG. protocols.io dx.doi.org/10.17504/protocols.io.uyjexun

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.uyjexun>

Authors: Alexandra Ehl, David Frommholz, Nadine Stefanczyk

Summary: Purification Guide for the Isolation of Antibodies with ChroPlate Filtration Plates by DALEX Biotech. Easy and quick high throughput antibody purification from various sources and species. Each well of the ChroPlate has a binding capacity of > 1 mg (tested with human polyclonal Ig, binding varies between species and clones). The proprietary resin does not shrink or swell in aqueous buffers. High pressure stability. pH stability short term 2 - 8, long term 3 - 8. Excellent thermal stability up to 15 minutes at 80 °C in aqueous buffers at neutral pH. Can be dried for long term storage (80 °C for > 2 h).

Affiliations: DALEX Biotech, DALEX Biotech, DALEX Biotech

External URL: <https://dalex-biotech.com/>

Version: 1

Publication Date: 2018

Proper Citation: Alexandra Ehl, David Frommholz, Nadine Stefanczyk 2018. ChroPlate - ProteinG. protocols.io dx.doi.org/10.17504/protocols.io.uyjexun

Record Creation Time: 20210329T220531+0000

Record Last Update: 20210329T220953+0000

Data and Source Information

Source: