

Resource Summary Report

Generated by [dkNET](#) on Aug 19, 2026

Anti-Tubulin beta-III, clone TU-20, Alexa Fluor 488 Conjugated

RRID:AB_1977541

Type: Antibody

Proper Citation

Millipore Cat# CBL412X, RRID:AB_1977541

Antibody Information

URL: http://antibodyregistry.org/AB_1977541

Proper Citation: Millipore Cat# CBL412X, RRID:AB_1977541

Target Antigen: Tubulin beta-III clone TU-20 Alexa Fluor 488

Host Organism: mouse

Clonality: monoclonal

Comments: seller recommendations: IgG1; IgG1 Immunofluorescence; Immunohistochemistry; Western Blot; Immunoprecipitation; Immunocytochemistry; IC, IF, WB, IH, IP

Antibody Name: Anti-Tubulin beta-III, clone TU-20, Alexa Fluor 488 Conjugated

Description: This monoclonal targets Tubulin beta-III clone TU-20 Alexa Fluor 488

Target Organism: b, porcine, h, m, r, pm, po

Antibody ID: AB_1977541

Vendor: Millipore

Catalog Number: CBL412X

Record Creation Time: 20250820T051958+0000

Record Last Update: 20250820T212011+0000

Ratings and Alerts

No rating or validation information has been found for Anti-Tubulin beta-III, clone TU-20, Alexa Fluor 488 Conjugated.

No alerts have been found for Anti-Tubulin beta-III, clone TU-20, Alexa Fluor 488 Conjugated.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 9 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Yada Y, et al. (2021) Establishment of induced pluripotent stem cells from schizophrenia discordant fraternal twins. Stem cell research, 55, 102504.

Romano E, et al. (2021) Generation of a homozygous CIITA knockout iPS cell line using the CRISPR-Cas9 system. Stem cell research, 57, 102580.

Trionfini P, et al. (2020) Generation of two isogenic knockout PKD2 iPS cell lines, IRFMNi003-A-1 and IRFMNi003-A-2, using CRISPR/Cas9 technology. Stem cell research, 42, 101667.

Sasaki-Honda M, et al. (2020) Generation of a transgene-free iPSC line and genetically modified line from a facioscapulohumeral muscular dystrophy type 2 (FSHD2) patient with SMCHD1 p.Lys607Ter mutation. Stem cell research, 47, 101884.

Romano E, et al. (2020) Generation of PKD1 mono-allelic and bi-allelic knockout iPS cell lines using CRISPR-Cas9 system. Stem cell research, 47, 101881.

Benedetti V, et al. (2018) Engineered Kidney Tubules for Modeling Patient-Specific Diseases and Drug Discovery. EBioMedicine, 33, 253.

Ciampi O, et al. (2018) Generation of two isogenic iPS cell lines (IRFMNi002-A and IRFMNi002-B) from a patient affected by Focal Segmental Glomerulosclerosis carrying a heterozygous c.565G>A mutation in PAX2 gene. Stem cell research, 33, 175.

Tan GW, et al. (2017) Induced pluripotent stem cells derived from an autosomal dominant lateral temporal epilepsy (ADLTE) patient carrying S473L mutation in leucine-rich glioma inactivated 1 (LGI1). Stem cell research, 24, 12.

Sekine SI, et al. (2017) Induced pluripotent stem cells derived from a patient with familial idiopathic basal ganglia calcification (IBGC) caused by a mutation in SLC20A2 gene. Stem cell research, 24, 40.