

Resource Summary Report

Generated by [dkNET](#) on Aug 20, 2026

Anti-Myc Tag, clone 4A6

RRID:AB_568800

Type: Antibody

Proper Citation

Millipore Cat# 05-724MG, RRID:AB_568800

Antibody Information

URL: http://antibodyregistry.org/AB_568800

Proper Citation: Millipore Cat# 05-724MG, RRID:AB_568800

Target Antigen: Myc Tag clone 4A6

Host Organism: mouse

Clonality: monoclonal

Comments: seller recommendations: IgG1; IgG1 Immunoprecipitation; Immunofluorescence; Immunocytochemistry; Western Blot; ChIP; ChIP, IF, IP, WB, IC

Antibody Name: Anti-Myc Tag, clone 4A6

Description: This monoclonal targets Myc Tag clone 4A6

Target Organism: h

Antibody ID: AB_568800

Vendor: Millipore

Catalog Number: 05-724MG

Record Creation Time: 20250819T224326+0000

Record Last Update: 20250820T030648+0000

Ratings and Alerts

No rating or validation information has been found for Anti-Myc Tag, clone 4A6.

No alerts have been found for Anti-Myc Tag, clone 4A6.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 9 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Nunez-Vazquez R, et al. (2025) The histone variant H3.14 is an early player in the abiotic stress response of Arabidopsis. *Developmental cell*.

Tahmasebinia F, et al. (2025) The 40S ribosomal subunit recycling complex modulates mitochondrial dynamics and endoplasmic reticulum - mitochondria tethering at mitochondrial fission/fusion hotspots. *Nature communications*, 16(1), 1021.

Osaka J, et al. (2024) Complex formation of immunoglobulin superfamily molecules Side-IV and Beat-IIb regulates synaptic specificity. *Cell reports*, 43(2), 113798.

Hernandez JC, et al. (2023) LIN28 and histone H3K4 methylase induce TLR4 to generate tumor-initiating stem-like cells. *iScience*, 26(3), 106254.

Shirasaki T, et al. (2022) Nonlytic Quasi-Enveloped Hepatovirus Release Is Facilitated by pX Protein Interaction with the E3 Ubiquitin Ligase ITCH. *Journal of virology*, 96(21), e0119522.

Harbauer AB, et al. (2022) Neuronal mitochondria transport Pink1 mRNA via synaptojanin 2 to support local mitophagy. *Neuron*, 110(9), 1516.

Kohashi K, et al. (2021) Sequential oncogenic mutations influence cell competition. *Current biology : CB*, 31(18), 3984.

Lee K, et al. (2021) Arabidopsis ATXR2 represses de novo shoot organogenesis in the transition from callus to shoot formation. *Cell reports*, 37(6), 109980.

Lee YC, et al. (2018) Znf179 E3 ligase-mediated TDP-43 polyubiquitination is involved in TDP-43- ubiquitinated inclusions (UBI) (+)-related neurodegenerative pathology. *Journal of biomedical science*, 25(1), 76.