

Resource Summary Report

Generated by [dkNET](#) on Aug 7, 2026

APC/Cyanine7 anti-human CD8

RRID:AB_2044006

Type: Antibody

Proper Citation

BioLegend Cat# 344714, RRID:AB_2044006

Antibody Information

URL: http://antibodyregistry.org/AB_2044006

Proper Citation: BioLegend Cat# 344714, RRID:AB_2044006

Target Antigen: CD8

Host Organism: mouse

Clonality: monoclonal

Comments: Applications: FC

Antibody Name: APC/Cyanine7 anti-human CD8

Description: This monoclonal targets CD8

Target Organism: cynomolgus, rhesus, human

Clone ID: Clone SK1

Antibody ID: AB_2044006

Vendor: BioLegend

Catalog Number: 344714

Alternative Catalog Numbers: 344713

Record Creation Time: 20250820T011408+0000

Record Last Update: 20250820T102346+0000

Ratings and Alerts

- Used for flow cytometry assay by the Human Islet Research Network community. Contact(s): [Jing Chen](#), [Mark Atkinson](#), [Todd Brusko](#), [Clayton Mathews](#) - Human Islets Research Network <https://hirnetwork.org/>

No alerts have been found for APC/Cyanine7 anti-human CD8.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 52 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Zhang Y, et al. (2026) Human receptive endometrial assembloid for deciphering the implantation window. eLife, 12.

Sajani A, et al. (2026) Heterogeneity and plasticity of the naive CD4+ T cell compartment. Cell reports, 45(3), 116980.

Li YR, et al. (2026) Protocol to generate human stem cell-derived CD70-directed allogeneic CAR-NKT cells for treating renal cell carcinoma. STAR protocols, 7(1), 104340.

Ruan Z, et al. (2026) Single-cell profiling of immune reset in patients with refractory generalized myasthenia gravis receiving autologous CD19/BCMA CAR-T cell therapy. Med (New York, N.Y.), 7(4), 101026.

Fan AE, et al. (2025) STAT5 Activation Enhances Adoptive Therapy Combined with Peptide Vaccination by Preventing PD-1 Inhibition. Molecular cancer therapeutics, 24(3), 419.

Schiantarelli J, et al. (2025) Genomic mediators of acquired resistance to immunotherapy in metastatic melanoma. Cancer cell, 43(2), 308.

Li YR, et al. (2025) Protocol to profile tumor and microenvironment from ovarian cancer patient samples and evaluate cell-based therapy using in vitro killing assays. STAR protocols, 6(2), 103742.

Webb L, et al. (2025) Membrane-bound IL-15 co-expression powers a potent and persistent CD70-targeted TRuC T-cell therapy. Frontiers in immunology, 16, 1609658.

Lyu Z, et al. (2025) Protocol for assessing pharmacokinetics and pharmacodynamics of human CAR-NKT cells in humanized mouse models using bioluminescence imaging. STAR protocols, 6(3), 103957.

Li B, et al. (2025) Serotonin transporter inhibits antitumor immunity through regulating the intratumoral serotonin axis. Cell, 188(14), 3823.

Li YR, et al. (2025) Overcoming ovarian cancer resistance and evasion to CAR-T cell therapy by harnessing allogeneic CAR-NKT cells. *Med (New York, N.Y.)*, 100804.

Thomé J, et al. (2025) C-reactive protein induced T cell activation is an indirect monocyte-dependent mechanism involving the CD80/CD28 pathway. *Frontiers in immunology*, 16, 1622865.

Gujar R, et al. (2025) Developing a therapeutic elastase that stimulates anti-tumor immunity by selectively killing cancer cells. *Cell reports. Medicine*, 6(11), 102446.

Kondo T, et al. (2025) Engineering TCR-controlled fuzzy logic into CAR T cells enhances therapeutic specificity. *Cell*, 188(9), 2372.

Li YR, et al. (2025) Generating allogeneic CAR-NKT cells for off-the-shelf cancer immunotherapy with genetically engineered HSP cells and feeder-free differentiation culture. *Nature protocols*, 20(5), 1352.

Song H, et al. (2025) Targeting tumor monocyte-intrinsic PD-L1 by rewiring STING signaling and enhancing STING agonist therapy. *Cancer cell*, 43(3), 503.

Fang Y, et al. (2025) Protocol to generate allojection-resistant universal CAR-NKT cells and evaluate their efficacy, mechanism of action, safety, and immunogenicity. *STAR protocols*, 6(2), 103810.

Semitekolou M, et al. (2025) Blood immunomap for prediction of responses to anti-PD-1 immunotherapy in metastatic non-small cell lung cancer. *iScience*, 28(9), 112804.

Weber J, et al. (2025) ROR2-specific CAR T cells are effective against hematologic and solid tumors and well tolerated in mice. *Cell reports. Medicine*, 6(10), 102400.

Kulkarni AJ, et al. (2025) L-Phenylalanine is a metabolic checkpoint of human Th2 cells. *Cell reports. Medicine*, 6(12), 102466.