

Resource Summary Report

Generated by [dkNET](#) on Aug 19, 2026

Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488

RRID:AB_2556606

Type: Antibody

Proper Citation

Thermo Fisher Scientific Cat# SA5-10026, RRID:AB_2556606

Antibody Information

URL: http://antibodyregistry.org/AB_2556606

Proper Citation: Thermo Fisher Scientific Cat# SA5-10026, RRID:AB_2556606

Target Antigen: Rat IgG (H+L)

Host Organism: donkey

Clonality: polyclonal secondary

Comments: Applications: ICC/IF (1:50-1:500), IHC (1:50-1:500), IP (Assay-dependent), WB (1:5,000-1:20,000), Flow (1:50-1:200)

Antibody Name: Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488

Description: This polyclonal secondary targets Rat IgG (H+L)

Target Organism: rat

Antibody ID: AB_2556606

Vendor: Thermo Fisher Scientific

Catalog Number: SA5-10026

Record Creation Time: 20250820T035327+0000

Record Last Update: 20250820T172723+0000

Ratings and Alerts

No rating or validation information has been found for Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488.

No alerts have been found for Donkey anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, DyLight™ 488.

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 7 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Nahiyera M, et al. (2025) Sca-1 expression depicts pro-inflammatory murine neutrophils under steady-state and pathological conditions. Cell reports, 44(7), 115995.

Puente-Ruiz SC, et al. (2025) B cell-derived nociceptin/orphanin FQ contributes to impaired glucose tolerance and insulin resistance in obesity. iScience, 28(7), 112819.

Zhang L, et al. (2024) Regulation of muscle hypertrophy through granulin: Relayed communication among mesenchymal progenitors, macrophages, and satellite cells. Cell reports, 43(4), 114052.

Sun F, et al. (2024) AdipoRon promotes amyloid-? clearance through enhancing autophagy via nuclear GAPDH-induced sirtuin 1 activation in Alzheimer's disease. British journal of pharmacology, 181(17), 3039.

Kaneshige A, et al. (2022) Relayed signaling between mesenchymal progenitors and muscle stem cells ensures adaptive stem cell response to increased mechanical load. Cell stem cell, 29(2), 265.

Winter C, et al. (2018) Chrono-pharmacological Targeting of the CCL2-CCR2 Axis Ameliorates Atherosclerosis. Cell metabolism, 28(1), 175.

Gabunia K, et al. (2016) IL-19 Halts Progression of Atherosclerotic Plaque, Polarizes, and Increases Cholesterol Uptake and Efflux in Macrophages. The American journal of pathology, 186(5), 1361.