

Resource Summary Report

Generated by [dkNET](#) on Sep 4, 2026

TSA™ Kit #15, with HRP—Goat Anti-Rabbit IgG and Alexa Fluor™ 594 Tyramide

RRID:AB_2716806

Type: Antibody

Proper Citation

Thermo Fisher Scientific Cat# T20925, RRID:AB_2716806

Antibody Information

URL: http://antibodyregistry.org/AB_2716806

Proper Citation: Thermo Fisher Scientific Cat# T20925, RRID:AB_2716806

Target Antigen: Rabbit

Host Organism: goat

Clonality: polyclonal

Comments: Discontinued; This kit contains one antibody: HRP—goat anti-rabbit IgG. Tyramide signal amplification is an enzyme-mediated detection method that uses the catalytic activity of horseradish peroxidase (HRP) to generate high-density labeling of a target protein or nucleic acid sequence in situ. The TSA Kit #15 includes HRP—goat anti-rabbit IgG and our bright, red-fluorescent Alexa Fluor 594 tyramide (which is spectrally similar to Texas Red dye).

Antibody Name: TSA™ Kit #15, with HRP—Goat Anti-Rabbit IgG and Alexa Fluor™ 594 Tyramide

Description: This polyclonal targets Rabbit

Target Organism: rabbit

Defining Citation: [PMID:16319317](#)

Antibody ID: AB_2716806

Vendor: Thermo Fisher Scientific

Catalog Number: T20925

Record Creation Time: 20231110T033803+0000

Record Last Update: 20260829T055218+0000

Ratings and Alerts

No rating or validation information has been found for TSA™ Kit #15, with HRP—Goat Anti-Rabbit IgG and Alexa Fluor™ 594 Tyramide.

Warning: Discontinued at Thermo Fisher Scientific
Discontinued; This kit contains one antibody: HRP—goat anti-rabbit IgG.
Tyramide signal amplification is an enzyme-mediated detection method that uses the catalytic activity of horseradish peroxidase (HRP) to generate high-density labeling of a target protein or nucleic acid sequence in situ. The TSA Kit #15 includes HRP—goat anti-rabbit IgG and our bright, red-fluorescent Alexa Fluor 594 tyramide (which is spectrally similar to Texas Red dye).

Data and Source Information

Source: [Antibody Registry](#)

Usage and Citation Metrics

We found 4 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Maloney R, et al. (2024) Efficacy and specificity of melanopsin reporters for retinal ganglion cells. The Journal of comparative neurology, 532(2), e25591.

Quattrochi LE, et al. (2019) The M6 cell: A small-field bistratified photosensitive retinal ganglion cell. The Journal of comparative neurology, 527(1), 297.

Itoga CA, et al. (2019) New viral-genetic mapping uncovers an enrichment of corticotropin-releasing hormone-expressing neuronal inputs to the nucleus accumbens from stress-related brain regions. The Journal of comparative neurology, 527(15), 2474.

Stabio ME, et al. (2018) The M5 Cell: A Color-Opponent Intrinsically Photosensitive Retinal Ganglion Cell. Neuron, 97(1), 150.