

Resource Summary Report

Generated by [dkNET](#) on Jul 30, 2026

Agilent: 6560 Ion Mobility LC/Q-TOF

RRID:SCR_019431

Type: Tool

Proper Citation

Agilent: 6560 Ion Mobility LC/Q-TOF (RRID:SCR_019431)

Resource Information

URL: <https://www.agilent.com/en/products/liquid-chromatography-mass-spectrometry-lc-ms/lc-ms-instruments/quadrupole-time-of-flight-lc-ms/6560-ion-mobility-q-tof-lc-ms>

Proper Citation: Agilent: 6560 Ion Mobility LC/Q-TOF (RRID:SCR_019431)

Description: Ion Mobility LC/Q-TOF delivers unrivaled separation and selectivity by combining chromatography, ion mobility, and mass spectrometry. Quadrupole Time-of-Flight liquid chromatography mass spectrometer that combines chromatography, ion mobility, and mass spectrometry. The system has electrodynamic ion funnel technology and an attractive low field drift tube design.

Synonyms: Agilent 6560 Ion Mobility Quadrupole Time-of-Flight Liquid Chromatography Mass Spectrometer

Resource Type: instrument resource

Keywords: Agilent, Quadrupole Time of Flight LC/MS, Instrument Equipment, USEdit,

Availability: Commercially available

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR_019431.pdf

Resource Name: Agilent: 6560 Ion Mobility LC/Q-TOF

Resource ID: SCR_019431

Alternate IDs: SCR_019434, Model_Number_6560

Record Creation Time: 20220129T080345+0000

Record Last Update: 20260728T043557+0000

Ratings and Alerts

No rating or validation information has been found for Agilent: 6560 Ion Mobility LC/Q-TOF.

No alerts have been found for Agilent: 6560 Ion Mobility LC/Q-TOF.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.