

# Resource Summary Report

Generated by [dkNET](#) on Jul 30, 2026

## Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator

RRID:SCR\_020859

Type: Tool

### Proper Citation

Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator (RRID:SCR\_020859)

### Resource Information

**URL:** <https://www.thermofisher.com/order/catalog/product/1567501?SID=srch-srp-1567501#/1567501?SID=srch-srp-1567501>

**Proper Citation:** Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator (RRID:SCR\_020859)

**Description:** Refrigerated/heated bath circulator that delivers high heating and cooling capacities for rapid heat-up and cool-down times and has electrical requirements of 200/230 V 60 Hz 200 V 50 Hz

**Resource Type:** instrument resource

**Keywords:** ThermoFisher, Refrigerated/Heated Bath Circulator, Instrument, Resource Equipment, Equipment, USEDit,

**Availability:** Commercially available

**Resource Name:** Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator

**Resource ID:** SCR\_020859

**Alternate IDs:** Model\_Number\_Ultra-low Temp G50 PC200 208-230 VOLT 60HZ

**Record Creation Time:** 20220129T080352+0000

**Record Last Update:** 20260725T190936+0000

### Ratings and Alerts

No rating or validation information has been found for Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator.

No alerts have been found for Thermo Fisher: Ultra-low Temp G50 PC200 208-230 VOLT 60HZ; Refrigerated/Heated Bath Circulator.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.