

Resource Summary Report

Generated by [dkNET](#) on Jul 29, 2026

Bertin Technologies: Precellys 24 tissue homogenizer

RRID:SCR_022979

Type: Tool

Proper Citation

Bertin Technologies: Precellys 24 tissue homogenizer (RRID:SCR_022979)

Resource Information

URL: <https://www.bertin-technologies.com/product/sample-preparation-homogenizers/precellys24-tissue-homogenizer/>

Proper Citation: Bertin Technologies: Precellys 24 tissue homogenizer (RRID:SCR_022979)

Description: The Precellys 24 homogenizer with automatic tube locking system with the ability to process up to 24 tubes simultaneously., THIS RESOURCE IS NO LONGER IN SERVICE. Documented on September 16,2025.

Synonyms: Bertin Precellys 24 tissue homogenizer

Resource Type: instrument resource

Keywords: USEDit, Instrument, Equipment, Bertin Technologies, homogenizer, automatic tube locking system,

Availability: THIS RESOURCE IS NO LONGER IN SERVICE

Resource Name: Bertin Technologies: Precellys 24 tissue homogenizer

Resource ID: SCR_022979

Alternate IDs: Model_Number_Precellys_24

Record Creation Time: 20221118T050202+0000

Record Last Update: 20260725T191026+0000

Ratings and Alerts

No rating or validation information has been found for Bertin Technologies: Precellys 24 tissue homogenizer.

No alerts have been found for Bertin Technologies: Precellys 24 tissue homogenizer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 3 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Tisi A, et al. (2025) Alterations of endocannabinoid signaling and microglia reactivity in the retinas of AD-like mice precede the onset of hippocampal β -amyloid plaques. Journal of neurochemistry, 169(2), e16256.

Oparija-Rogenmozere L, et al. (2025) Time-restricted feeding leads to sex- and organ-specific responses in the murine digestive system. PloS one, 20(9), e0332295.

Zapata-Pavas LE, et al. (2025) Protocol for pull-down capture of ubiquitinated and NEDDylated liver proteins from transgenic biotinylated ubiquitin and NEDD8 mice. STAR protocols, 6(2), 103697.