

Resource Summary Report

Generated by [dkNET](#) on Aug 6, 2026

Covaris: CP02 Automated Dry Pulverizer

RRID:SCR_026012

Type: Tool

Proper Citation

Covaris: CP02 Automated Dry Pulverizer (RRID:SCR_026012)

Resource Information

URL: <https://www.covaris.com/products-services/instruments/cryoprep>

Proper Citation: Covaris: CP02 Automated Dry Pulverizer (RRID:SCR_026012)

Description: Covaris cryoPREP Dry Pulverizer enables dry, non-contact cryopulverization of tissue samples in a closed system without the risk of carryover or cross-contamination. The cryoPREP system provides a simple and robust tissue processing method for low and high mass inputs. The specially engineered tissueTUBE is used to process samples and extract biomolecules from a diverse set of tissue samples ranging in mass from 15 mg up to 2 g. Both manual and automated versions of the cryoPREP Dry Pulverizers are available to suit your needs.

Synonyms: Covaris cryoPREP Dry Pulverizers

Resource Type: instrument resource

Keywords: automated, dry, non-contact cryopulverization, tissue samples, closed system,

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/main/PDF/SCR_026012.pdf

Resource Name: Covaris: CP02 Automated Dry Pulverizer

Resource ID: SCR_026012

Record Creation Time: 20241101T053303+0000

Record Last Update: 20260801T191341+0000

Ratings and Alerts

No rating or validation information has been found for Covaris: CP02 Automated Dry Pulverizer.

No alerts have been found for Covaris: CP02 Automated Dry Pulverizer.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We have not found any literature mentions for this resource.