

Resource Summary Report

Generated by [dkNET](#) on Aug 17, 2026

FITMASTER

RRID:SCR_016233

Type: Tool

Proper Citation

FITMASTER (RRID:SCR_016233)

Resource Information

URL: http://www.heka.com/downloads/downloads_main.html#down_fitmaster

Proper Citation: FITMASTER (RRID:SCR_016233)

Description: Software for analysis and fitting routines of electrophysiological data. Analysis can be performed on the levels of Sweeps/Traces and Series.

Synonyms: HEKA FITMASTER

Resource Type: software resource, software application, data processing software, data analysis software

Keywords: electrophysiology, analysis, fitting, statistic, polynomial, exponential, gaussian, boltzmann, formalism, voltage, current

Availability: Free, Available for download, Tutorial available, Demo available

Resource Name: FITMASTER

Resource ID: SCR_016233

License URLs: http://www.heka.com/about/about_main.html#terms

Record Creation Time: 20220129T080329+0000

Record Last Update: 20260817T040548+0000

Ratings and Alerts

No rating or validation information has been found for FITMASTER.

No alerts have been found for FITMASTER.

Data and Source Information

Usage and Citation Metrics

We found 101 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Falck A, et al. (2026) Potassium-chloride cotransporter 2 activity dampens induced ictal-like activity in neocortical slices containing the seizure propagation zone of temporal lobe epilepsy patients. *Epilepsia*, 67(1), 499.

Lesnak JB, et al. (2026) An improved preservation method for human dorsal root ganglion neurons enables wider access to human molecular pain neuroscience. *Cell reports methods*, 6(5), 101412.

Cordeiro S, et al. (2026) Natural xanthenes as γ -Mangostin induce vasorelaxation involving key gating residues in the S6 domain of BK channels. *eLife*, 14.

Ghovanloo MR, et al. (2025) SCN3A-related neurodevelopmental disorder: Clinical case reports and biophysical characterization. *Channels (Austin, Tex.)*, 19(1), 2580175.

Centonze E, et al. (2025) Case Report: Functional investigation of the γ ENaC G532S mutation presenting as mild PHA-1B3. *Frontiers in medicine*, 12, 1605057.

Sesti V, et al. (2025) Membrane-targeted push-pull azobenzenes for the optical modulation of membrane potential. *Light, science & applications*, 14(1), 8.

Fan X, et al. (2025) Exosomal mir-126-3p derived from endothelial cells induces ion channel dysfunction by targeting RGS3 signaling in cardiomyocytes: a novel mechanism in Takotsubo cardiomyopathy. *Stem cell research & therapy*, 16(1), 36.

Stojanovic S, et al. (2025) Rebound Bursting Selectively Enables Fast Dynamics in Dopamine Midbrain Neurons Projecting to the Dorsolateral Striatum. *The Journal of neuroscience : the official journal of the Society for Neuroscience*, 45(44).

Pampols-Perez M, et al. (2025) Mechanosensitive PIEZO2 channels shape coronary artery development. *Nature cardiovascular research*, 4(7), 921.

Moubarak E, et al. (2025) Postnatal Development of Dendritic Morphology and Action Potential Shape in Rat Substantia Nigra Dopaminergic Neurons. *eNeuro*, 12(4).

Zhang J, et al. (2025) The transcription factor ZFP64 promotes activity-dependent synapse elimination during postnatal cerebellar development. *iScience*, 28(6), 112746.

Hartveit E, et al. (2025) HCN channels in rod bipolar cells of rat retina: subcellular localization, kinetic properties and functional dynamics. *The Journal of physiology*.

Skrabak D, et al. (2023) Slack K⁺ channels limit kainic acid-induced seizure severity in mice by modulating neuronal excitability and firing. *Communications biology*, 6(1), 1029.

Kapell H, et al. (2023) Neuron-oligodendrocyte potassium shuttling at nodes of Ranvier protects against inflammatory demyelination. *The Journal of clinical investigation*, 133(7).

Hussein RA, et al. (2023) Selenomethionine mis-incorporation and redox-dependent voltage-gated sodium channel gain of function. *Journal of neurochemistry*, 167(2), 262.

König C, et al. (2023) Prostaglandin EP3 receptor activation is antinociceptive in sensory neurons via PI3K γ , AMPK and GRK2. *British journal of pharmacology*, 180(4), 441.

Verkest C, et al. (2022) Intrinsically disordered intracellular domains control key features of the mechanically-gated ion channel PIEZO2. *Nature communications*, 13(1), 1365.

El Hamdaoui Y, et al. (2022) Analysis of hyperforin (St. John's wort) action at TRPC6 channel leads to the development of a new class of antidepressant drugs. *Molecular psychiatry*, 27(12), 5070.

Lemmens V, et al. (2022) Hetero-pentamerization determines mobility and conductance of Glycine receptor α 3 splice variants. *Cellular and molecular life sciences : CMLS*, 79(11), 540.

Dai XQ, et al. (2022) Heterogenous impairment of β cell function in type 2 diabetes is linked to cell maturation state. *Cell metabolism*, 34(2), 256.