

# Resource Summary Report

Generated by [dkNET](#) on Sep 12, 2026

## Agilent: 8800 Triple Quadrupole ICP-MS

RRID:SCR\_019457

Type: Tool

---

### Proper Citation

Agilent: 8800 Triple Quadrupole ICP-MS (RRID:SCR\_019457)

---

### Resource Information

**URL:** [https://www.agilent.com/cs/library/applications/appcompendium\\_icp-qqq-5991-2802en-us-agilent.pdf](https://www.agilent.com/cs/library/applications/appcompendium_icp-qqq-5991-2802en-us-agilent.pdf)

**Proper Citation:** Agilent: 8800 Triple Quadrupole ICP-MS (RRID:SCR\_019457)

**Description:** Triple Quadrupole ICP-MS enables MS MS operation, providing control of reaction processes in the ORS3 collision and reaction cell.

**Resource Type:** instrument resource

**Keywords:** Agilent, ICP-MS, Instrument Equipment, USEdit,

**Availability:** Restricted

**Specification URL:** [https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR\\_019457.pdf](https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR_019457.pdf)

**Resource Name:** Agilent: 8800 Triple Quadrupole ICP-MS

**Resource ID:** SCR\_019457

**Alternate IDs:** Model\_Number\_Agilent\_8800\_ICP-QQQ

**Alternate URLs:** <https://hpst.cz/sites/default/files/oldfiles/agilent-8800-icp-qqq-application-handbook.pdf>

**Old URLs:** <https://www.agilent.com/en/products/icp-ms/icp-ms-systems/8800-triple-quadrupole-icp-ms>

**Record Creation Time:** 20220129T080345+0000

**Record Last Update:** 20260912T195915+0000

---

### Ratings and Alerts

No rating or validation information has been found for Agilent: 8800 Triple Quadrupole ICP-MS.

No alerts have been found for Agilent: 8800 Triple Quadrupole ICP-MS.

---

## Data and Source Information

**Source:** [SciCrunch Registry](#)

---

## Usage and Citation Metrics

We have not found any literature mentions for this resource.