

Resource Summary Report

Generated by [dkNET](#) on Sep 15, 2026

Perkin Elmer: Opera Phenix Plus High Content Screening System

RRID:SCR_021100

Type: Tool

Proper Citation

Perkin Elmer: Opera Phenix Plus High Content Screening System (RRID:SCR_021100)

Resource Information

URL: <https://resources.revvity.com/pdfs/tch-opera-phenix-plus-improved-3d-imaging.pdf>

Proper Citation: Perkin Elmer: Opera Phenix Plus High Content Screening System (RRID:SCR_021100)

Description: Opera Phenix Plus system for high throughput high content assays, phenotypic screening, assays using complex disease models, such as live cells, primary cells and microtissues, and fast response assays like calcium flux or cardiomyocyte beating. Optical design lets you generate information through confocal imaging and at higher throughput through simultaneous acquisition. Delivers speed without compromising sensitivity. Uses Harmony High Content Imaging and Analysis Software to quantify complex cellular phenotypes. Harmony software is designed for PerkinElmer high content screening systems.

Synonyms: Opera Phenix Plus

Resource Type: instrument resource

Keywords: PerkinElmer, confocal imaging, simultaneous acquisition, high throughput high content assays, phenotypic screening, Harmony software, high content screening systems, instrument, equipment, USEDit

Availability: Commercially available

Specification URL: https://raw.githubusercontent.com/SciCrunch/RRID-Instruments/refs/heads/main/PDF/SCR_021100.pdf

Resource Name: Perkin Elmer: Opera Phenix Plus High Content Screening System

Resource ID: SCR_021100

Alternate IDs: Model_Number_Opera_Phenix_Plus

Old URLs: <https://www.perkinelmer.com/Product/opera-phenix-plus-system-hh14001000>

Record Creation Time: 20220129T080353+0000

Record Last Update: 20260912T195951+0000

Ratings and Alerts

No rating or validation information has been found for Perkin Elmer: Opera Phenix Plus High Content Screening System.

No alerts have been found for Perkin Elmer: Opera Phenix Plus High Content Screening System.

Data and Source Information

Source: [SciCrunch Registry](#)

Usage and Citation Metrics

We found 14 mentions in open access literature.

Listed below are recent publications. The full list is available at [dkNET](#) .

Hough SH, et al. (2026) Loss of CTLH component MAEA impairs DNA repair and replication and leads to developmental delay. EMBO molecular medicine, 18(2), 492.

Trautmann RK, et al. (2026) High-Throughput 3D Glioblastoma Model in Glycosaminoglycan Hydrogels for Personalized Therapeutic Screening. Macromolecular bioscience, 26(1), e00394.

Fieblinger T, et al. (2026) Repurposing drugs for treating the neurobehavioral manifestations of PTEN hamartoma tumor syndrome. iScience, 29(4), 115428.

Kafka E, et al. (2026) TCA cycle rewiring underpins histone acetylation sourcing and cell-fate transitions during exit from naive pluripotency. Cell stem cell, 33(5), 820.

Li R, et al. (2026) Colony-stimulating factor-2 triggers NNMT-dependent myofibroblast activation in head and neck cancer. Oncogene, 45(16), 1425.

Shen F, et al. (2025) Human Kallikrein 2: A Novel Lineage-Specific Surface Target in Prostate Cancer. Clinical cancer research : an official journal of the American Association for Cancer Research, 31(21), 4543.

Weidling I, et al. (2025) hiPSC-neurons recapitulate the subtype-specific cell intrinsic nature of susceptibility to neurodegenerative disease-relevant aggregation. Acta neuropathologica communications, 13(1), 108.

Haladik B, et al. (2025) Image-based drug screening combined with molecular profiling identifies signatures and drivers of therapy resistance in pediatric AML. Cell reports.

Medicine, 6(9), 102304.

Lyu TT, et al. (2025) The sodium-glucose co-transporter 2 inhibitor, empagliflozin, attenuates pulmonary vascular remodelling by inhibiting the phosphorylation of PDGF receptor-?. British journal of pharmacology.

Kortleve D, et al. (2024) TCR-Engineered T Cells Directed against Ropporin-1 Constitute a Safe and Effective Treatment for Triple-Negative Breast Cancer. Cancer discovery, 14(12), 2450.

Badja C, et al. (2024) Insights from multi-omic modeling of neurodegeneration in xeroderma pigmentosum using an induced pluripotent stem cell system. Cell reports, 43(6), 114243.

Lee H, et al. (2023) In vitro characterization on the role of APOE polymorphism in human hippocampal neurogenesis. Hippocampus, 33(4), 322.

Berryer MH, et al. (2023) Robust induction of functional astrocytes using NGN2 expression in human pluripotent stem cells. iScience, 26(7), 106995.

Kodani Y, et al. (2022) EpCAM Is a Surface Marker for Enriching Anterior Pituitary Cells From Human Hypothalamic-Pituitary Organoids. Frontiers in endocrinology, 13, 941166.