

Resource Summary Report

Generated by [dkNET](#) on Sep 16, 2026

BiomekFXp Robot Minipreps (RoboPreps)

DOI:10.17504/protocols.io.3zqgp5w

Type: Protocol

Proper Citation

James Angstman 2019. BiomekFXp Robot Minipreps (RoboPreps). protocols.io
dx.doi.org/10.17504/protocols.io.3zqgp5w

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.3zqgp5w>

Authors: James Angstman

Group: The Center for Genome Editing and Recording

Summary: RoboPreps are great for high-throughput cloning and subcloning, especially if you're going to be doing Maxipreps afterward anyway. RoboPreps must be done in 96 well format as of now, though future generations may one day break through this technological barrier.

Affiliations: Harvard University

Version: 1

Publication Date: 2019

Proper Citation: James Angstman 2019. BiomekFXp Robot Minipreps (RoboPreps). protocols.io dx.doi.org/10.17504/protocols.io.3zqgp5w

Record Creation Time: 20210329T220517+0000

Record Last Update: 20210329T220835+0000

Data and Source Information

Source: