

Resource Summary Report

Generated by [dkNET](#) on Sep 18, 2026

TXTL phage production

DOI:10.17504/protocols.io.77thrnn

Type: Protocol

Proper Citation

Marijn Ceelen 2019. TXTL phage production. protocols.io
dx.doi.org/10.17504/protocols.io.77thrnn

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.77thrnn>

Authors: Marijn Ceelen

Group: iGEM Wageningen 2019

Summary: This protocol explains how to use the Abor BiosciencesSigma 70 TXTL mix to produce bacteriophages. This protocol works for bacteriophage T7. It has been tested on phage lambda. However, it seems that the current protocol does not work for phage lambda.

Affiliations: Wageningen University

Version: 1

Publication Date: 2019

Proper Citation: Marijn Ceelen 2019. TXTL phage production. protocols.io
dx.doi.org/10.17504/protocols.io.77thrnn

Record Creation Time: 20210329T220523+0000

Record Last Update: 20210329T220901+0000

Data and Source Information

Source: