

Resource Summary Report

Generated by [dkNET](#) on Sep 9, 2026

Phage lambda genome concatemerization

DOI:10.17504/protocols.io.7k3hkyn

Type: Protocol

Proper Citation

Marijn Ceelen 2019. Phage lambda genome concatemerization. protocols.io
[dx.doi.org/10.17504/protocols.io.7k3hkyn](https://doi.org/10.17504/protocols.io.7k3hkyn)

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.7k3hkyn>

Authors: Marijn Ceelen

Group: iGEM Wageningen 2019

Summary: This protocol is used to create concatemers of phage genomes. This is done by ligating phage lambda genomes using the overhangs at the cos site with T4 ligase.

Affiliations: Wageningen University

Version: 1

Publication Date: 2019

Proper Citation: Marijn Ceelen 2019. Phage lambda genome concatemerization. protocols.io [dx.doi.org/10.17504/protocols.io.7k3hkyn](https://doi.org/10.17504/protocols.io.7k3hkyn)

Record Creation Time: 20210329T220533+0000

Record Last Update: 20210329T221013+0000

Data and Source Information

Source: