

Resource Summary Report

Generated by [dkNET](#) on Sep 17, 2026

Phosphate Buffered Saline (PBS)

DOI:10.17504/protocols.io.bezgjf3w

Type: Protocol

Proper Citation

Allen Institute for Brain Science 2020. Phosphate Buffered Saline (PBS). protocols.io
dx.doi.org/10.17504/protocols.io.bezgjf3w

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.bezgjf3w>

Authors: Allen Institute for Brain Science

Group: BICCN, Allen Institute for Brain Science

Summary: This protocol is used to prepare Phosphate Buffered Saline (PBS). PBS (-) and PBS (+) which includes Tween are used as biologically neutral wash solutions in many processes including In Situ Hybridization (ISH), Immunohistochemistry (IHC) and Fixation. Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Affiliations: Allen Institute

Version: 1

Publication Date: 2020

Proper Citation: Allen Institute for Brain Science 2020. Phosphate Buffered Saline (PBS). protocols.io dx.doi.org/10.17504/protocols.io.bezgjf3w

Record Creation Time: 20210329T220524+0000

Record Last Update: 20210329T220903+0000

Data and Source Information

Source: