

Resource Summary Report

Generated by [dkNET](#) on Sep 6, 2026

Preparation of acute midbrain slices for patch-clamp recordings

DOI:10.17504/protocols.io.bhfzj3p6

Type: Protocol

Proper Citation

Oriol Pavón Arocas, Tiago Branco 2021. Preparation of acute midbrain slices for patch-clamp recordings. protocols.io dx.doi.org/10.17504/protocols.io.bhfzj3p6

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.bhfzj3p6>

Authors: Oriol Pavón Arocas, Tiago Branco

Summary: Here we describe two different approaches to prepare acute midbrain slices from young adult mice for slice electrophysiology experiments. We also provide general advice for troubleshooting common issues leading to suboptimal slice quality, as well as some suggestions to ensure good maintenance of a patch-clamp rig.

Affiliations: Sainsbury Wellcome Centre for Neural Circuits and Behaviour, University College London, Howland Street, London, W1T 4JG, UK, Sainsbury Wellcome Centre for Neural Circuits and Behaviour, University College London, Howland Street, London, W1T 4JG, UK

Version: 1

Publication Date: 2021

Proper Citation: Oriol Pavón Arocas, Tiago Branco 2021. Preparation of acute midbrain slices for patch-clamp recordings. protocols.io dx.doi.org/10.17504/protocols.io.bhfzj3p6

Record Creation Time: 20210329T220528+0000

Record Last Update: 20210329T220933+0000

Data and Source Information

Source: