

Resource Summary Report

Generated by [dkNET](#) on Sep 16, 2026

Chimeric Protein-LAG and Peptostreptococcal protein L sandwich ELISA

DOI:10.17504/protocols.io.btbaniie

Type: Protocol

Proper Citation

Angel Justiz-Vaillant 2021. Chimeric Protein-LAG and Peptostreptococcal protein L sandwich ELISA. protocols.io dx.doi.org/10.17504/protocols.io.btbaniie

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.btbaniie>

Authors: Angel Justiz-Vaillant

Group: Carbon

Summary: This ELISA is used to study the interactions between protein-LAG (PLAG) and Peptostreptococcal protein-L (SpL) with different immunoglobulin preparations from mammalian and avian species.

Affiliations: University of the West Indies St. Augustine

Version: 1

Publication Date: 2021

Proper Citation: Angel Justiz-Vaillant 2021. Chimeric Protein-LAG and Peptostreptococcal protein L sandwich ELISA. protocols.io dx.doi.org/10.17504/protocols.io.btbaniie

Record Creation Time: 20210329T220516+0000

Record Last Update: 20210329T220831+0000

Data and Source Information

Source: