

Resource Summary Report

Generated by [dkNET](#) on Sep 16, 2026

Resolving gel buffer (4X)

DOI:10.17504/protocols.io.eq3bdyn

Type: Protocol

Proper Citation

DAVID DUNIGAN AND IRINA AGARKOVA 2016. Resolving gel buffer (4X). protocols.io
dx.doi.org/10.17504/protocols.io.eq3bdyn

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.eq3bdyn>

Authors: DAVID DUNIGAN AND IRINA AGARKOVA

Group: VERVE Net

Summary: For use in Polyacrylamide Gel System For Electrophoresis Of Proteins.

Version: 1

Publication Date: 2016

Proper Citation: DAVID DUNIGAN AND IRINA AGARKOVA 2016. Resolving gel buffer (4X). protocols.io dx.doi.org/10.17504/protocols.io.eq3bdyn

Record Creation Time: 20210329T220516+0000

Record Last Update: 20210329T220830+0000

Data and Source Information

Source: