

Resource Summary Report

Generated by [dkNET](#) on Sep 8, 2026

Fluorescent Focus Unit Assay using LICOR Imaging System

DOI:10.17504/protocols.io.k8fcztn

Type: Protocol

Proper Citation

Bernardo Mainou 2017. Fluorescent Focus Unit Assay using LICOR Imaging System. protocols.io dx.doi.org/10.17504/protocols.io.k8fcztn

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.k8fcztn>

Authors: Bernardo Mainou

Summary: Indirect immunofluorescence infectivity assay for reovirus using LICOR Imaging System

Affiliations: Emory University

External URL: <https://doi.org/10.1371/journal.ppat.1006768>

Version: 2

Publication Date: 2017

Proper Citation: Bernardo Mainou 2017. Fluorescent Focus Unit Assay using LICOR Imaging System. protocols.io dx.doi.org/10.17504/protocols.io.k8fcztn

Record Creation Time: 20210329T220530+0000

Record Last Update: 20210329T220946+0000

Data and Source Information

Source: