

Resource Summary Report

Generated by [dkNET](#) on Sep 20, 2026

Chromatin immunoprecipitation

DOI:10.17504/protocols.io.knkcvcw

Type: Protocol

Proper Citation

Gabriel Sanchez, Rosemary Kiernan 2018. Chromatin immunoprecipitation. protocols.io dx.doi.org/10.17504/protocols.io.knkcvcw

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.knkcvcw>

Authors: Gabriel Sanchez, Rosemary Kiernan

Summary: This protocol can be used for chromatin immunoprecipitation of RNAPII and associated factors, as well as histones. The settings are given for HeLa cells and should be adapted for other cell types.??

Associated Publications: Contreras X, Salifou K, Sanchez G, Helsmoortel M, Beyne E, Bluy L, Pelletier S, Rousset E, Rouquier S, Kiernan R, Nuclear RNA surveillance complexes silence HIV-1 transcription. PLoS Pathogens 14(3). doi: [10.1371/journal.ppat.1006950](https://doi.org/10.1371/journal.ppat.1006950)

Affiliations: Institut de Genetique Humaine, UMR9002, Montpellier, France, Institut de Genetique Humaine, UMR9002, Montpellier, France

External URL: <https://doi.org/10.1371/journal.ppat.1006950>

Version: 2

Publication Date: 2018

Proper Citation: Gabriel Sanchez, Rosemary Kiernan 2018. Chromatin immunoprecipitation. protocols.io dx.doi.org/10.17504/protocols.io.knkcvcw

Record Creation Time: 20210329T220527+0000

Record Last Update: 20210329T220927+0000

Data and Source Information

Source: