

Resource Summary Report

Generated by [dkNET](#) on Sep 6, 2026

Colony PCR (DreamTaq) with Synechocystis sp. PCC 6803 and gel electrophoresis

DOI:10.17504/protocols.io.mk5c4y6

Type: Protocol

Proper Citation

Maximilian Dietsch, Miriam D 2018. Colony PCR (DreamTaq) with Synechocystis sp. PCC 6803 and gel electrophoresis. protocols.io dx.doi.org/10.17504/protocols.io.mk5c4y6

Protocol Information

URL: <https://dx.doi.org/10.17504/protocols.io.mk5c4y6>

Authors: Maximilian Dietsch, Miriam D

Group: Axmann Lab, CyanoWorld, iGEM Duesseldorf 2018

Summary: In this protocol, the extraction of template DNA out of Synechocystis sp. PCC 6803 samples, the PCR of the template and the application on a agarosegel are described. The protocol was handed over by Maximilian Dietsch MSc.

Affiliations: Institute for Synthetic Microbiology, Heinrich-Heine Universität Düsseldorf

Version: 1

Publication Date: 2018

Proper Citation: Maximilian Dietsch, Miriam D 2018. Colony PCR (DreamTaq) with Synechocystis sp. PCC 6803 and gel electrophoresis. protocols.io dx.doi.org/10.17504/protocols.io.mk5c4y6

Record Creation Time: 20210329T220528+0000

Record Last Update: 20210329T220936+0000

Data and Source Information

Source: